

Product Datasheet

Aminoallyl-dUTP (6-JOE) (orb64170)

Description	Aminoallyl-dUTP conjugated to 6-JOE.
Form/Appearance	solution in water; pH: 7.5 ± 0.5
Concentration	1.0 mM-1.1 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ _{exc} 520 nm, λ _{em} 548 nm, ε 75.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₃₅ H ₃₂ Cl ₂ N ₃ O ₂₂ P ₃
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1010.47 g/mol (free acid); Detected MW: 1010.01 g/mol (free acid)
SMILES	C(/C=C/c1cn(c(=O)[nH]c1=O)[C@H]1C[C@@H]([C@H](O1)COP(=O)(O)OP(=O)(O)OP(=O)(O)O)NC(=O)c1cc2C3(c4cc(c(O)c(Cl)c4O)c4c(Cl)c(O)c(cc34)OC)OC)OC(=O)c2cc1
Expiration Date	12 months from date of receipt.