

Product Datasheet

CSDE1 antibody (orb641791)

Description	Rabbit polyclonal antibody to CSDE1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Recombinant protein of human CSDE1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 102 kDa
Uniprot ID	O75534
Dilution Range	WB1:500 - 1:2000IP1:20 - 1:50
Expiration Date	12 months from date of receipt.