

Product Datasheet

HIF2a antibody (orb641997)

Description	Rabbit polyclonal antibody to HIF2a
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant Protein of human HIF2a.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 120 kDa
Uniprot ID	Q99814
Dilution Range	WB 1:500 - 1:2000, IF 1:50 - 1:100
Expiration Date	12 months from date of receipt.