

Product Datasheet

PRICKLE2 antibody (orb642064)

Description	Rabbit polyclonal antibody to PRICKLE2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	KLH conjugated synthetic peptide derived from human Prickle 2:401-500/844
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 96 kDa
Uniprot ID	Q7Z3G6
Dilution Range	WB 1:500 - 1:2000, IF 1:50 - 1:100
Expiration Date	12 months from date of receipt.