

Product Datasheet

THOC1 antibody (orb642079)

Description	Rabbit polyclonal antibody to THOC1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant protein of human THOC1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 92 kDa
Uniprot ID	Q96FV9
Dilution Range	WB 1:500 - 1:2000 , IHC-P: -P/, IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.