

Product Datasheet

PPP4R1 antibody (orb642104)

Description	Rabbit polyclonal antibody to PPP4R1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant protein of human PPP4R1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 130 kDa
Uniprot ID	Q8TF05
Dilution Range	WB 1:500 - 1:2000 , IHC-P: -P 1:50 - 1:200IP 1:50 - 1:100
Expiration Date	12 months from date of receipt.