

Product Datasheet

Beclin 1 antibody (orb642253)

Description	Rabbit polyclonal antibody to Beclin 1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fragment corresponding to a region within amino acids 141-451 of Human Beclin 1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 55 kDa
Uniprot ID	Q14457
Dilution Range	WB: 1:2000-1:5000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.