

Product Datasheet

ATP1A1 (phospho-S23) antibody (orb642305)

Description	Rabbit polyclonal antibody to ATP1A1 (phospho-S23)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human ATP1A1 around the phosphorylation site of Serine23.
Preservatives	1 mg/ml in Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 100 kDa
Uniprot ID	P05023
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.