

Product Datasheet

TUBA1A (phospho-Y272) antibody (orb642311)

Description	Rabbit polyclonal antibody to TUBA1A (phospho-Y272)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human TUBA1A around the phosphorylation site of Tyrosine272.
Preservatives	1 mg/ml in Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 55 kDa
Uniprot ID	Q6PEY2 , Q9BQE3 , P68363 , Q71U36 , P68366
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.