

Product Datasheet

NMDA Epsilon 2 (phospho-Y1336) antibody (orb642320)

Description	Rabbit polyclonal antibody to NMDA Epsilon 2 (phospho-Y1336)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human NMDA ϵ 2 around the phosphorylation site of Tyrosine1336.
Preservatives	1 mg/ml in Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 190 kDa
Uniprot ID	Q13224
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.