

Product Datasheet

HA-tag (4G3) antibody (orb642367)

Description	Mouse monoclonal antibody to HA-tag (4G3)
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide HA-tag.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB:1:1000-1:5000 , IF:1:1000 IP:1:100
Expiration Date	12 months from date of receipt.