

Product Datasheet

Smad3 (phospho-T179) antibody (orb642405)

Description	Rabbit polyclonal antibody to Smad3 (phospho-T179)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human Smad3 around the phosphorylation site of Threonine 179.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 48, 55 kDa
Uniprot ID	P84022
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.