

Product Datasheet

GFP-tag (3A10) antibody (orb642448)

Description	Mouse monoclonal antibody to GFP-tag (3A10)
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Synthetic peptide GFP-tag.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB:1:5000-1:10000IP: 1:100-500
Expiration Date	12 months from date of receipt.