

Product Datasheet

VEGF-A (L98) antibody (orb642460)

Description	Rabbit polyclonal antibody to VEGF-A (L98)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 70-120 of Human VEGF206.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 27 kDa (VEGF206); 45 kDa (L-VEGF206); ~ 15-44 kDa (Others Isoform)
Uniprot ID	P15692
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.