

Product Datasheet

E2F-1 (D426) antibody (orb642826)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to E2F-1 (D426) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC |
| Immunogen | Synthetic peptide, corresponding to amino acids 400-450 of Human E2F-1. |
| Preservatives | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2 |
| Concentration | 1mg/ml |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 47 kDa |
| Uniprot ID | Q01094 |
| Dilution Range | IHC-P: 1:50-1:200 |
| Expiration Date | 12 months from date of receipt. |