

Product Datasheet

IRS-1 (V317) antibody (orb642867)

Description	Rabbit polyclonal antibody to IRS-1 (V317)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Synthetic peptide, corresponding to amino acids 300-350 of Human IRS-1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~132, 180 kDa
Uniprot ID	P35568
Dilution Range	IHC-P: 1:50-1:200 , IF: 1:50-1:200
Expiration Date	12 months from date of receipt.