

Product Datasheet

Rsk-1 (T353) antibody (orb642870)

Description	Rabbit polyclonal antibody to Rsk-1 (T353)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Synthetic peptide, corresponding to Human Rsk-1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 85 kDa
Uniprot ID	Q15418
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.