

Product Datasheet

Tubulin beta (4D9) antibody (orb642897)

Description	Mouse monoclonal antibody to Tubulin beta (4D9)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 402-455 of Human Tubulin β .
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 55 kDa
Uniprot ID	P68371, P04350, Q13885, Q9BUF5, Q13509, P07437
Dilution Range	WB: 1:5000-20000, IHC-P: 1:100-500, IF: 1:100-500
Expiration Date	12 months from date of receipt.