

Product Datasheet

MAX protein (phospho-Tyr123) antibody (orb6430)

Description	MAX protein (phospho-Tyr123) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MAX protein around the phosphorylation site of Tyr123 (SL(p-Y)TN)
Target	phospho-MAX protein (Tyr123)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18 kDa
Uniprot ID	P61244
Expiration Date	12 months from date of receipt.