

Product Datasheet

EphA2/3/4 (T587) antibody (orb643019)

Description	Rabbit polyclonal antibody to EphA2/3/4 (T587)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic peptide, corresponding to amino acids 550-600 of Human EphA2.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 112 kDa
Uniprot ID	P29317, P29320, P54764
Dilution Range	WB: 1:500-1:1000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.