

Product Datasheet

Actinin- alpha 1/2/3/4 (W48) antibody (orb643402)

Description	Rabbit polyclonal antibody to Actinin-alpha 1/2/3/4 (W48)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 30-60 of Human Actinin- α 1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 103 kDa
Uniprot ID	O43707 , Q08043 , P12814 , P35609
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.