

Product Datasheet

AR- alpha 1A (Q365) antibody (orb643403)

Description	Rabbit polyclonal antibody to AR- alpha 1A (Q365)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 350-400 of Human AR- α 1A.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 51 kDa
Uniprot ID	P35348
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.