

Product Datasheet

TAF I p48 (R318) antibody (orb643481)

Description	Rabbit polyclonal antibody to TAF I p48 (R318)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 280-336 of Human TAF I p48.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 53 kDa
Uniprot ID	Q15573
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.