

Product Datasheet

IRAK-M (L519) antibody (orb643737)

Description	Rabbit polyclonal antibody to IRAK-M (L519)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Synthetic peptide, corresponding to amino acids 500-540 of Human IRAK-M.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 68 kDa
Uniprot ID	Q9Y616
Dilution Range	IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.