

Product Datasheet

MARK1/2/3/4 (N211) antibody (orb643849)

Description	Rabbit polyclonal antibody to MARK1/2/3/4 (N211)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Synthetic peptide, corresponding to amino acids 180-230 of Human MARK1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 80 to 92 kDa
Uniprot ID	P27448 , Q9P0L2 , Q7KZI7 , Q96L34
Dilution Range	IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.