

Product Datasheet

CHP2 (M136) antibody (orb643969)

Description	Rabbit polyclonal antibody to CHP2 (M136)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide, corresponding to amino acids 122-182 of Human MARK1/2/3/4.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 23 kDa
Uniprot ID	O43745
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.