

Product Datasheet

ERK1/2 (D196) antibody (orb644157)

Description	Rabbit polyclonal antibody to ERK1/2 (D196)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic peptide, corresponding to amino acids 164-209 of Human ERK1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 42,44 kDa
Uniprot ID	P27361, P28482
Dilution Range	WB: 1:500-1:1000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.