

Product Datasheet

ZIP-kinase (K261) antibody (orb644232)

Description	Rabbit polyclonal antibody to ZIP-kinase (K261)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide, corresponding to amino acids 231-280 of Human ZIP-kinase.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 52 kDa
Uniprot ID	O43293
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.