

Product Datasheet

LAT (D156) antibody (orb644255)

Description	Rabbit polyclonal antibody to LAT (D156)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 131-180 of Human LAT.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 36, 38 kDa
Uniprot ID	O43561
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.