

Product Datasheet

Elk1 (phospho-S383) antibody (orb644401)

Description	Rabbit polyclonal antibody to Elk1 (phospho-S383)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Synthetic phosphopeptide derived from human Elk-1 around the phosphorylation site of Serine 383.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 45 kDa
Uniprot ID	P19419
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.