

Product Datasheet

Integrin beta 3 (phospho-Y773) antibody (orb644429)

Description	Rabbit polyclonal antibody to Integrin beta 3 (phospho-Y773)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic phosphopeptide derived from human Integrin β 3 around the phosphorylation site of Tyrosine 773.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 87, 110 kDa
Uniprot ID	P05106
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.