

Product Datasheet

JAK1 (phospho-Y1022) antibody (orb644432)

Description	Rabbit polyclonal antibody to JAK1 (phospho-Y1022)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic phosphopeptide derived from human JAK1 around the phosphorylation site of Tyrosine 1022.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 125 kDa
Uniprot ID	P23458
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.