

Product Datasheet

MEF-2 (phospho-T312) antibody (orb644441)

Description	Rabbit polyclonal antibody to MEF-2 (phospho-T312)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	Synthetic phosphopeptide derived from human MEF-2 around the phosphorylation site of Threonine 312.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 55 kDa
Uniprot ID	Q02078
Dilution Range	IHC-P: 1:50-1:200IP: 1:10-1:100
Expiration Date	12 months from date of receipt.