

Product Datasheet

Stat3 (phospho-Y705) antibody (orb644483)

Description	Rabbit polyclonal antibody to Stat3 (phospho-Y705)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Synthetic phosphopeptide derived from human Stat3 around the phosphorylation site of Tyrosine 705.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 88 kDa
Uniprot ID	P40763
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200IP: 1:50-1:200
Expiration Date	12 months from date of receipt.