

Product Datasheet

Tau (phospho-T181) antibody (orb644494)

Description	Rabbit polyclonal antibody to Tau (phospho-T181)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic phosphopeptide derived from human isoforms Tau-4 around the phosphorylation site of Threonine 181.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	50 to 80 kDa
Uniprot ID	P10636
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.