

Product Datasheet

Dynamin-1 (phospho-S774) antibody (orb644518)

Description	Rabbit polyclonal antibody to Dynamin-1 (phospho-S774)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human Dynamin-1 around the phosphorylation site of Serine 774.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 98 kDa
Uniprot ID	Q05193
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.