

Product Datasheet

IKK Gamma (phospho-S31) antibody (orb644525)

Description	Rabbit polyclonal antibody to IKK Gamma (phospho-S31)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human IKK γ around the phosphorylation site of Serine 31.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 48 kDa
Uniprot ID	Q9Y6K9
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.