

Product Datasheet

Na⁺/K⁺-ATPase alpha 1 (phospho-S16) antibody (orb644539)

Description	Rabbit polyclonal antibody to Na ⁺ /K ⁺ -ATPase alpha 1 (phospho-S16)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic phosphopeptide derived from human Na ⁺ /K ⁺ -ATPase α1 around the phosphorylation site of Serine 16.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 113 kDa
Uniprot ID	P05023
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.