

Product Datasheet

JNK1/2/3 (phospho-T183/Y185) antibody (orb644571)

Description	Rabbit polyclonal antibody to JNK1/2/3 (phospho-T183/Y185)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Synthetic phosphopeptide derived from human JNK1 around the phosphorylation site of Threonine 183 and Tyrosine 185.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 46,54 kDa
Uniprot ID	P53779, P45983, P45984
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.