

Product Datasheet

Catenin- Delta 1 (phospho-Y228) antibody (orb644592)

Description	Rabbit polyclonal antibody to Catenin-Delta 1 (phospho-Y228)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human Catenin- δ 1 around the phosphorylation site of Tyrosine 228
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 95, 100 kDa
Uniprot ID	O60716
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.