

Product Datasheet

Ron (phospho-S1394) antibody (orb644593)

Description	Rabbit polyclonal antibody to Ron (phospho-S1394)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human Ron around the phosphorylation site of Serine 1394
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 145, 150 kDa
Uniprot ID	Q04912
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.