

Product Datasheet

14-3-3 Theta /Tau (phospho-S232) antibody (orb644614)

Description	Rabbit polyclonal antibody to 14-3-3 Theta /Tau (phospho-S232)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human 14-3-3 θ/τ around the phosphorylation site of Serine 232
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -22°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 28 kDa
Uniprot ID	P27348
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.