

## Product Datasheet

### NIPA (phospho-S354) antibody (orb644627)

|                            |                                                                                                |
|----------------------------|------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to NIPA (phospho-S354)                                              |
| <b>Species/Host</b>        | Rabbit                                                                                         |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                                   |
| <b>Tested Applications</b> | WB                                                                                             |
| <b>Immunogen</b>           | Synthetic phosphopeptide derived from human NIPA around the phosphorylation site of Serine 354 |
| <b>Preservatives</b>       | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2                        |
| <b>Concentration</b>       | 1mg/ml                                                                                         |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.       |
| <b>Note</b>                | For research use only                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                                     |
| <b>MW</b>                  | ~ 60-65 kDa                                                                                    |
| <b>Uniprot ID</b>          | <b>Q86WB0</b>                                                                                  |
| <b>Dilution Range</b>      | WB: 1:500-1:1000                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                |