

Product Datasheet

JNK1/2/3 (phospho-Y185) antibody (orb644702)

Description	Rabbit polyclonal antibody to JNK1/2/3 (phospho-Y185)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human JNK1 around the phosphorylation site of Tyrosine 185.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 46,54 kDa
Uniprot ID	P53779, P45983, P45984
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.