

Product Datasheet

HSF1 (phospho-S303) antibody (orb644723)

Description	Rabbit polyclonal antibody to HSF1 (phospho-S303)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Synthetic phosphopeptide derived from human HSF1 around the phosphorylation site of Serine 303.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 82 kDa
Uniprot ID	Q00613
Dilution Range	IHC-P: 1:50-1:200 , IF: 1:50-1:200
Expiration Date	12 months from date of receipt.