

Product Datasheet

Smad3 (phospho-S204) antibody (orb644765)

Description	Rabbit polyclonal antibody to Smad3 (phospho-S204)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	Synthetic phosphopeptide derived from human Smad3 around the phosphorylation site of Serine 204.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 48, 55 kDa
Uniprot ID	P84022
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200ICC: 1:50-1:200IP: 1:50-1:200
Expiration Date	12 months from date of receipt.