

Product Datasheet

VEGFR2 (phospho-Y1175) antibody (orb644792)

Description	Rabbit polyclonal antibody to VEGFR2 (phospho-Y1175)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 150, 210 kDa
Uniprot ID	P35968
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.