

Product Datasheet

SLC6A1 antibody (orb644944)

Description	Rabbit polyclonal antibody to SLC6A1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 431-476 of Human SLC6A1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 67 kDa
Uniprot ID	P30531
Dilution Range	WB: 1:500-1:1000ICC: 1:50-1:200
Expiration Date	12 months from date of receipt.