

Product Datasheet

COL26A1 antibody (orb645355)

Description	Rabbit polyclonal antibody to COL26A1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to residues in Human COL26A1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 45, 55 kDa
Uniprot ID	Q96A83
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.