

Product Datasheet

CCL26 antibody (orb645511)

Description	Rabbit polyclonal antibody to CCL26
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to residues in HumanCCL26.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -23°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 11 kDa
Uniprot ID	Q9Y258
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:203
Expiration Date	12 months from date of receipt.